


Cleaning Your ITE (In The Ear) Hearing Aid with Wax Tube

How to keep your hearing aid working

Wax blocking the sound outlet of the hearing aid is a very common reason for aids to stop working. Sound simply cannot get out of the hearing aid and into your ear.

A simple cleaning procedure every time you change your battery will prevent this.



1. Locate the wax tube, which collects wax (see above).
2. Hold the aid with the wax tube pointed downward.
3. Using a wax pick, which is a small plastic rod with a wire hook, gently remove all of the wax that is in the wax tube. This will allow your aid to work properly.